


Simon Reed^{1,2}, Patrick van Eijk^{1,2} & Felix Dobbs²

Background |

Advances in genome editing are facilitating the development of novel cell and gene therapies. These innovations enable precise modifications to existing genes and the creation of synthetic genes, particularly beneficial for generating immune-compatible cells for allogeneic therapy, providing readily available treatment options. Despite the potential of gene editing in disease treatment, ensuring its safe and effective application demands the development of new tools. Risks associated with genome editing include the potential for malignant transformation of cultured cells, attributed to mutations acquired during cell culture, as well as the induction of genomic instability, especially in cells with pre-existing DNA repair defects. Moreover, the tools used for gene editing can inadvertently cause DNA breaks at unintended sites in the genome, leading to off-target effects. Such off-target editing poses risks such as the activation of proto-oncogenes or disruption of tumor suppressors, potentially driving carcinogenesis. Additionally, off-target mutagenesis may generate neoantigens, triggering autoimmune reactions or other forms of cellular dysfunction. Identifying and addressing off-target gene editing is crucial in mitigating these risks.

Figure 1 | Gene-editing is a two-step process



Step 1 introduces some form of damage in the DNA targeted to a specific location (e.g. a DSB induced by CRISPR-Cas9). Step 2, gene editing occurs as a result of endogenous DNA repair mechanism in the cell, typically utilising one of two major DNA repair pathways (Non-Homologous End Joining (NHEJ) and Homology-Direct Repair (HDR)). The outcome of the repair reaction results in a mutation at the target-site, modifying the gene of interest.

INDUCE-seq™ |


More precise methods are needed for testing off-target gene-editing during all phases of therapeutic development, including treatment follow-up. At present, standardised assays to assess the safety of gene editing-based therapies are lacking¹. INDUCE-seq was developed to address this². INDUCE-seq is a scalable platform technology for mapping and characterizing DNA breaks. It leverages a novel PCR-free methodology for *in situ* break capture and sequencing by NGS, revealing the breaks induced by any nuclease-based genome editing system with high precision. INDUCE-seq is the first unbiased cell-based solution that is free from PCR induced biases that distort measurements, has broad compatibility with a wide range of therapeutically relevant cells, and applicable to any nuclease-based gene editing system. INDUCE-seq provides data-driven and actionable insights to accelerate research & development, pre-clinical and clinical stages gene editing programs.

INDUCE-SEQ PERFORMANCE METRICS

INDUCE-seq's ability to differentiate between induced vs endogenous recurrent breaks was calculated using a standard induced break dataset.

Metrics	Description	Value
Sensitivity	How well we can call editing-induced breaks in a background of endogenous ones	~94%
Specificity	Our ability to call endogenous breaks from induced ones at low recurrence	~85%
Precision	Our ability to accurately call induced breaks at high recurrency	~96%
FNR	Occur at a recurrency at which mutations occur at background levels	~6.5%
FPR	What fraction of our off-target list of recurrent breaks are false positives	~15%

Figure 2 | Detecting editing-induced breaks using INDUCE-seq




In situ break labelling in fixed and permeabilised cells is performed by ligating a full-length, chemically modified P5 sequencing adapter to end-prepared DSBs. Genomic DNA is then extracted, fragmented, end-prepared and ligated using a chemically modified half-functional P7 adapter. Resulting DNA libraries contain a mixture of functional DSB-labelled fragments (P5:P7) and non-functional genomic DNA fragments (P7:P7). Subsequent sequencing of INDUCE-seq libraries enriches for DSB-labelled fragments and eliminates all other non-functional DNA. As the INDUCE-seq library preparation is PCR-free, each sequencing read obtained is equivalent to a single labelled DSB-end from a cell².

HESI Cross-Validation Study |

Here we describe a pilot study to evaluate the method involving a consortium of partners under the umbrella of HESI CT-TRACS committee³. CRISPR-Cas9-based gene editing of five well-studied genetic targets was conducted by two independent industry partners (AZ & Novartis), using two different cell types. On- and off-target gene editing was assessed by measuring breaks in the genome using INDUCE-seq. The genetic changes at these locations were subsequently measured using error-corrected sequencing using Duplex-seq⁴ which allows for sensitive detection of mutations.

Figure 3 | Identifying on & off-target editing-induced breaks with INDUCE-seq



Top Panel | Using the unique digital read-out of the INDUCE-seq DSB detection method, we developed an off-target nomination process that combines simple break recurrency with guide sequence homology to generate a superset of potential off-targets. **Middle Panel |** A hierarchical break report is generated enabling a data-driven analysis of off-target likelihood, as shown in three classes (i) common, (ii) unbiased & (iii) sequence-based. **Bottom Panel |** Genome browser view demonstrating break recurrency at HBD gene-edit locus.

Figure 4 | The experimental plan undertaken by the CT-TRACS consortium





Figure 5 | Mutation analysis at on- and off-target sites by DuplexSeq

AAVS1 Off-target Report: Mutation & Break Data Summary

Off-target name	chr	start	end	Break Score		Editing Efficiency					
				GM24385	MCF10A	24h	48h	Control	24h	48h	Control
01_AAVS1_common	19	55,627,121	55,627,127	673	673	59.7%	72.56%	0.09%	88.76%	83.45%	0.04%
02_AAVS1_common	6	36,765,477	36,765,483	40	156	12.5%	7.8%	0.00%	27%	22%	0.00%
03_AAVS1_common	13	106,612,923	106,612,929	1	2	3.1%	5.89%	0.00%	2.93%	2.90%	0.00%
04_AAVS1_common	12	107,486,278	107,486,282	16	24	1.7%	2.98%	0.00%	2.91%	3.74%	0.00%
05_AAVS1_common	5	169,911,417	169,911,421	17	23	2.5%	4.74%	0.00%	5.34%	6.10%	0.01%
06_AAVS1_common	11	61,111,251	61,111,254	16	8	0.4%	0.71%	0.00%	1.63%	3.29%	0.00%
07_AAVS1_common	10	132,580,828	132,580,833	12	3	0.1%	0.13%	0.00%	0.29%	0.42%	0.00%
08_AAVS1_common	8	144,885,132	144,885,138	10	4	0.8%	0.13%	0.00%	1.17%	0.86%	0.00%
09_AAVS1_common	8	22,633,566	22,633,569	6	4	1.23%	1.60%	0.00%	2.24%	1.81%	0.00%
10_AAVS1_common	15	67,137,278	67,137,283	8	1	0.01%	0.00%	0.00%	0.02%	0.04%	0.00%
11_AAVS1_common	11	62,297,767	62,297,771	8	0	0.01%	0.00%	0.00%	0.04%	0.03%	0.00%
12_AAVS1_common	10	121,188,694	121,188,697	4	3	0.09%	0.16%	0.00%	0.26%	0.29%	0.00%
13_AAVS1_common	6	53,200,090	53,200,095	6	0	0.09%	0.19%	0.01%	0.03%	0.03%	0.00%
14_AAVS1_common	15	69,476,691	69,476,696	4	1	0.02%	0.03%	0.01%	0.02%	0.07%	0.00%
15_AAVS1_common	21	42,892,941	42,892,964	6	5	0.28%	0.54%	0.00%	0.69%	1.38%	0.00%
16_AAVS1_common	21	29,716,114	29,716,137	6	3	0.81%	0.01%	0.00%	0.97%	0.99%	0.00%
17_AAVS1_common	22	22,479,063	22,479,086	7	2	0.01%	0.00%	0.00%	0.03%	0.27%	0.00%
18_AAVS1_common	10	116,704,562	116,704,585	4	3	0.01%	0.02%	0.00%	0.19%	0.17%	0.00%
19_AAVS1_common	1	49,443,700	49,443,723	7	6	0.00%	0.01%	0.00%	0.02%	0.02%	0.00%
20_AAVS1_common	13	100,289,739	100,289,762	4	2	0.00%	0.00%	0.00%	0.01%	0.01%	0.00%
21_AAVS1_common	16	31,484,902	31,484,920	3	3	0.00%	0.01%	0.00%	0.01%	0.01%	0.00%
22_AAVS1_common	6	24,124,802	24,124,815	3	3	0.00%	0.01%	0.00%	0.01%	0.00%	0.00%
23_AAVS1_common	15	89,175,079	89,175,102	2	3	0.02%	0.05%	0.04%	0.08%	0.09%	1.76%
24_AAVS1_common	2	205,396,126	205,396,149	4	1	0.21%	0.38%	0.00%	0.86%	1.24%	0.00%
25_AAVS1_common	12	108,581,675	108,581,698	4	0	0.00%	0.00%	0.00%	0.07%	0.10%	0.00%
26_AAVS1_common	15	49,827,892	49,827,915	3	1	0.01%	0.02%	0.00%	0.04%	0.07%	0.00%
27_AAVS1_common	17	77,121,114	77,121,137	4	6	0.00%	0.00%	0.00%	0.00%	0.01%	0.00%
28_AAVS1_common	12	49,461,976	49,461,983	5	0	0.01%	0.00%	0.00%	0.04%	0.01%	0.00%
29_AAVS1_common	1	162,106,706	162,106,709	11	0	0.01%	0.00%	0.00%	0.01%	0.00%	0.00%
30_AAVS1_common	2	70,475,841	70,475,841	11	0	0.01%	0.00%	0.00%	0.00%	0.01%	0.00%
31_AAVS1_common	2	159,047,922	159,047,923	11	0	0.00%	0.00%	0.00%	0.00%	0.01%	0.00%
32_AAVS1_common	5	68,630,719	68,630,726	0	0	0.00%	0.00%	0.00%	0.00%	0.01%	0.00%

Mutations at the high-confidence class of off-target are shown here, demonstrating a strong proportional relationship between break count and the accumulation of mutations. Importantly, mutation rates in the control sample indicate the background level of mutation.

Figure 6 | INDUCE-seq reproducibility



INDUCE-seq measures DSBs in an unbiased way. An inventory of the high-level break numbers are presented here. **A) Total Breaks:** Breaks measured by INDUCE-seq are reproducible between triplicates. GM24385 and MCF10A cells have different total amounts of DSBs that accumulate over time post gene-editing. **B) Number of Off-targets:** The number of sequence-based off-targets detected by INDUCE-seq is consistent between replicates. The majority of CRISPR-Cas9 activity is registered at 6 hours after transfection. Off-targets discovery is consistent between the two experiments. **C) On- & Off-target Break relationship:** Detailed analysis of the on- and off-target DSBs demonstrates a highly reproducible break detection across the 5 gene-editing experiments. All data shown plotted as the mean for each triplicate with error bars representing the standard deviation.

Conclusions |

- INDUCE-seq reproducibly detects on-target editing between the two laboratories and cell types.
- INDUCE-seq reproducibly identifies off-targets between laboratories.
- INDUCE-seq detects both common and unique off-targets in different cell types.
- DSB recurrency at both on- and off-targets correlates with mutational frequency.